



HydroMARK™


Breast Biopsy Site Marker

Being Different Can Make
All the Difference




H₂O

Made of exclusive hydrogel technology, the polyethylene glycol compounds attract water molecules creating the anechoic appearance under ultrasound.



12 months

The HydroMARK™ marker carrier, hydrogel, hydrates to provide up to 12 months of unequivocal ultrasound visibility.¹



MRI T2

The same hydrogel is uniquely visible under the T2 sequence where it presents as a bright white artifact.²

The HydroMARK™ marker is strikingly anechoic under ultrasound and uniquely distinguishable under T2 weighted MRI.^{1,2}

Easy. High Quality. Long Lasting. It's so different, you won't miss it.

Devicor Medical Europe GmbH

Carl-Zeiss-Straße 9 | 25451 Quickborn (Germany) **Phone:** +49 (0) 4106 978 9300 | **Customer Service:** 0800 500 98 98 | **Email:** info-europe@mammotome.com
Commercial Register: HRB 333205 | District court: Kiel | Tax number: 18/297/28192 | VAT-ID: DE202370516

1. Sakamoto, N., Fukuma, E., Tsunoda, Y. et al. Breast Cancer 25, 575–582 (2018). <https://doi.org/10.1007/s12282-018-0854-8>

2. Shah AD, Mehta AK, Talati N, Brem R, Margolies LR. Clin Imaging. 2019 May-Jun;55:196-212. doi: 10.1016/j.clinimag.2019.05.014. Epub 2019 Jun 14. PMID: 31204189.

Product may not be approved or available in your region. Please check with your local Mammotome representative.

© 2023 Devicor Medical Products, Inc. All rights reserved. Devicor, Mammotome and HydroMARK are registered trademarks of Devicor Medical Products, Inc. in the USA and optionally in other countries. Other logos, product and/or company names might be trademarks of their respective owners. MDM# 230232 Rev 3/23

Mammotome
mammotome.com