

MammoSTAR®

Biopsy Site Marker

A robust marker for sensitive application



Long-Term Ultrasound Visibility

Provides enhanced long-term ultrasound visibility to allow easier pre-surgery localization.¹




Axillary Lymph Node Application

Indicated for soft tissue including breast tissue and axillary lymph nodes.






All-Natural Composition

Offers a non-metallic option composed of all natural minerals, Carbon Coated Zirconium Oxide and lyophilized beta glucan gel (or simple sugar) in regular and petite.



Ceramic marker with
Pyrolytic Carbon Coating

Lyophilized Beta
Glucan Gel

Product Code	Gauge	Description	Qty Per Box
STAR1401	14G	3mm x 1mm barbell shape within beta glucan gel 	10
STAR1402	14G	5mm x 1mm tribell shape within beta glucan gel 	10
STAR1403	14G	3mm x 1mm barbell shape within petite beta glucan gel 	10

Devicor Medical Europe GmbH

Carl-Zeiss-Straße 9 | 25451 Quickborn (Germany) **Phone:** +49 (0) 4106 978 9300 | **Customer Service:** 0800 500 98 98 | **Email:** info-europe@mammotome.com
Commercial Register: HRB 333205 | District court: Kiel | Tax number: 18/297/28192 | VAT-ID: DE202370516

¹. Mammotome Internal Data

Long-term ultrasound visibility image courtesy of Dr. Kimberly C. Hutcherson, North Metropolitan Radiology Associates, LLP

Axillary lymph node images courtesy of Dr. Nicole Nelson, St. Rita's Medical Center

Product may not be approved or available in your region. Please check with your local Mammotome representative.

© 2023 Devicor Medical Products, Inc. All rights reserved. Devicor, Mammotome and MammoSTAR are registered trademarks of Devicor Medical Products, Inc., in the USA and optionally in other countries. Other logos, product and/or company names might be trademarks of their respective owners. MDM# 230233 Rev 3/23

Mammotome
mammotome.com